Select the degree of eye contact

Close eye contact

- Hold the diopter ring. Rotate the eyecup anticlockwise to release it.
- Pull out the eyecup.
- Hold the diopter ring. Rotate the eyecup clockwise to tighten it.

If you are a *spectacle wearer*, you need the whole of the available distance. Move the eyecups to the lowest position.

Set the magnification scale

You can set the magnification scale for various optics combinations:

- For the top range, 25× to 250×
- For the bottom range, 2.5× to 25×.

Example:
1.0× objective, 20× eyepiece

- Turn the magnification changer fully towards you (lowest magnification).
- Use the table on page 39 to determine the magnification range (16× to 160× in this example).

- Hook a finger into the back window of the magnification changer and displace the scale so that the number 16 is opposite to the index.

Use this procedure to set the scales on *both rotary knobs*.

Change the magnification

You can alter the magnification steplessly within the range 1:10 by means of the rotary knobs on each side of the stereomicroscope. (See magnification table on page 40.)

Engage the click stops

Turn the rotary knob towards the side-faced column until the magnification changer engages in *nine positions*.

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
Adjust the dioptic setting

Adjust the dioptic setting exactly for each eye.

If, before working with the instrument, each user carries out this procedure exactly, the image will be parfocal, i.e. it will remain sharp when the magnification is changed.

A Preparations

1. Adjust the illumination (pages 20–23).
2. Set the interpupillary distance and the viewing height (page 10).
3. Open the double-iris diaphragm if one is fitted.
4. Adjust the eyecups (page 11).
5. Set the zero diopter position on each eyepiece as shown in the diagram.

B Focus on the test object

6. Select the lowest magnification.
7. Position a flat test object beneath the objective.
8. Look into the eyepieces and use the coarse focusing knob to bring the image into focus.
9. Select the highest magnification.
10. Use the fine focusing knob to perfect the focusing.

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
C Adjust the dioptric setting

11. Select the *lowest magnification*.
12. Do not look into the eyepieces. Rotate the eyepieces anticlockwise (in the "+" direction) as far as the stop (+5 diopters).

13. Look into the eyepieces. Slowly rotate each eyepiece individually (in the "-" direction) until each eye sees the object sharply imaged.

14. Reset the eyepieces to "0". Repeat steps 9 to 13 twice more and make a note of the dioptric settings. Finally, set the median values on the eyepieces.

Inspect the parfocality

15. Select the *highest magnification*.
16. If necessary, refocus using the fine focusing drive.

17. Slowly move the magnification changer from the lowest magnification to the highest.

The image sharpness must remain constant.

Graticule (reticle) in eyepiece?

1. Rotate the dioptric ring of the eyepiece with graticule from "+" to "-" until you see the graticule sharply imaged.
2. Adjust the eyepiece without graticule as described above.

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Assembly:
Basic outfits with
incident- and
transmitted-light
stands

Incident-light stand

Secure the side-faced column:

- Introduce the three hollow screws into the baseplate *from beneath* and secure them to the end of the side-faced column.
- Fit the stage plate and the stage clips.

Transmitted-light stands

Secure the side-faced column:

- Pass the three hollow screws through the *adapter* and secure them to the end of the side-faced column.

- Use two hollow screws to secure the side-faced column to the *transmitted-light base*.

Further details can be found in the user manuals for the transmitted-light stands.

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M10 con estativo de diascopia EB

Optics carrier/ drive housing

Slide the optics carrier, together with the drive housing, *gently* down the side-faced column until you can use the focusing drive.

Objective

Screw on the selected objective *anticlockwise*.

If access to the object is inadequate, unscrew and remove the *adapter ring* used for fitting accessories to the 1.6× objective.

Binocular tube

- Release the clamping screw.
- Introduce the accessory tube into the dovetail ring and *displace it slightly from side to side* until the locating pin engages the groove.

Eyepieces

- Fit the graticule if required (page 33).
- Release the clamping screws of the eyepieces.
- Insert the eyepieces *fully* into the eyepiece tubes.
- Tighten the screws.

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M10 with incident-light stand

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M10 con estativo de episcopio

Swinging-arm and table-clamp stands

Tilt the side-faced column

- Release the clamping lever.
- Tilt the instrument.
- Tighten the clamping lever.

All clamping levers can be pulled out along their axes and freely repositioned.

Change the working distance and swing out laterally

- Release the clamping screw.
- Raise or lower the horizontal arm, or swing it out.
- Retighten the clamping screw.
- Use the safety ring.

Change the cantilever (overhang)

- Release the clamping screw.
- Displace the horizontal arm.
- Tighten the clamping screw.

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
M10 with swinging-arm stand

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Assembly

Use the clamp of the *table clamp stand* to secure it to tables between 20mm and 50mm thick.

Pass the 4 screws through the base of the *swinging-arm stand* from beneath and secure them to the end of the column.

- Secure the *safety ring* on the column with its screw head uppermost.

The safety ring, which ensures that the instrument does not slip down, must *always* be fitted beneath the horizontal arm.

- Lower the *horizontal arm* on to the safety ring so that the screw head engages the groove.

Tighten the clamping screw.

- Slacken the clamping lever (a).
- Unfold the *side-faced carrier rod* as shown in the diagram.
- Tighten the clamping lever.
- Secure the *inclined side-faced carrier rod* to the horizontal arm by means of the clamping lever (b).
- Unscrew and remove the safety plate.
- Slide the stereomicroscope, together with the drive housing, on to the side-faced column from beneath, until you can use the focusing drive.

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Universal stand

In this outfit, the optics carrier and the drive housing for 50mm diameter columns are together secured to the universal stand.

Assembly

- Introduce three hollow screws into the base-plate from beneath through the holes on the narrow side and secure them to the end of the column.

- Secure the safety ring to the column, screw head upwards.

The safety ring prevents the instrument from slipping down and must always be secured beneath the drive housing.

- Position the drive housing on the safety ring so that the screw head engages the groove.

- Release the two laterally-positioned hollow screws and extract the connecting plate from the drive housing.

- Position the optics carrier and the connecting plate as shown in the diagram.
- Introduce two hollow screws into the connecting holes as shown in the diagram and secure them finger-tight.

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




– Level up the connecting plate, using a spirit level or a set-square, and tighten the hollow screws.

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– Slide the optics carrier, together with the connecting plate, into the guide of the drive housing, and secure it with the two lateral hollow screws.

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Stage carriers with magnetic linkages can be attached to the baseplate.

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Fitting the cross stage
Introduce spacing pieces between the cross stage and the stage carrier, and use hollow screws to connect all of the components.

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M10 with universal stand and with 6in×6in cross stage

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**Incident illuminators:
6V/20W low-voltage lamp**

The colour temperature (3200K) is ideal for photography.

Adjust the lamp

- Engage the condenser lens at the click-stop position.
- Switch on the illumination (page 23).
- Hold the lamp housing and release the clamping screw.
- Align the lamp with the object, and tighten the clamping screw.
- To change the size of the spot, push the condenser lens forwards or backwards.

Filters

The filterholder with 32mm diameter diffusing filter can be clipped to the pegs along with three other filters.

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