## Possibilities for fitting

- on the cast foot
- on side-faced columns
- on the transmitted-light stand EB
- on the rear of the optics carrier (for the discussion tube)
- on the drive housing for 50mm diameter columns.

## Fitting

- Hold the lamp and release the clamping screw.
- Screw the thread on to one of the components listed above.
- Tighten the clamping screw.
- Connect the transformer (page 23).

## Changing the bulb

- Press the spring clips together and pull off the rear of the housing.
- Hold the new halogen bulb in a cloth and push it into position.
- Replace the rear of the lamp housing.

## Centre the bulb

Every time after a bulb change!

- Switch on the lamp at the transformer (page 23).
- Swing out the diffusing filter.
- Pull out the condenser lens completely.
- Hold the lamp and release the clamping screw.
- Direct the lamp vertically downwards on to a white paper so that two filament images are visible.
- Tighten the clamping screw.
- Focus the filament images as well as possible, using the slotted screw.
- Displace the filament images so that they are adjacent to one another and just in contact.

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**Incident illuminators:  
Coaxial incident illuminator**

With this illuminator you can study and photograph *flat highly-reflecting objects* such as wafers, integrated circuits, and polished metal sections.

The coaxial incident light housing fits like an additional tube *between the optics carrier and the binocular tube* or between the optics carrier and the phototube (page 24).

Please refer to the *detailed instructions* in the user manual M2 232.

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
M10 with coaxial incident light housing and with fibre-optic light guide

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## Regulating transformer

The 0–7V/50VA regulating transformer will accept the 6V/20W low-voltage lamp together with the 6V/20W coaxial incident illuminator.

Two lamps, with a maximum combined loading of 50W, can be connected to the sockets.

### Set for the voltage available

Displace the voltage selector on the underside:

- Position 115 for power supply 100V–120V
- Position 220 for power supply 200V–240V.

### Set the click stop at 6V

This prevents overloading of the lamps.

- Connect the power cable. Do not connect the lamp.
- Turn the regulating knob clockwise until it engages. If the indicator now shows 6V, the setting is correct.

If the indicator does not show 6V:

- Release the screw.
- Turn the regulating knob anticlockwise until the indicator is at zero.
- Connect the lamp required to one of the two sockets.
- Turn the regulating knob until the indicator is at 6V.
- Tighten the screw on the metal ring.

### Changing the fuse

Pull out the power cable.

- Remove two screws from the base.
- Lift off the housing.
- Fit a new 8A, 5×20 fuse (214477).
- Close the housing.

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




## Fit the accessory tubes

### Fitting tubes between optics carrier and binocular tube:

- Release the clamping screw.
- Introduce the accessory tube into the dovetail ring as already described for the binocular tube (page 15) and *displace it slightly from side to side* until the locating pin engages the groove.
- Tighten the clamping screw.
- Fit the binocular tube to the accessory tube in the same way.

### Filter slide housing

- Secure the housing to the optics carrier.
- Cut  $23\text{mm} \times 45\text{mm}$  strips from the gelatine filters and insert them into the filter slide.
- Push the filter slide into its housing.

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

M10 with trinocular tube and with WILD MPS46/52 with 35mm magazine

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### Phototubes

Phototubes serve for fitting the Wild Modular Photomicrographic Systems (MPS).

Please read the user manual for your photomicrographic system.

If you are not taking photographs, engage the observation beam path (VIS). When the photographic beam path (PHOT) is engaged, part of the light is directed into the camera.

### Trinocular tube

When you pull out the dark slide lever to the left, the dark slide stops light from entering or leaving the right eyepiece tube.

### Fitting an SLR camera

1. Insert the eyepiece tube into the phototube and secure it.
2. Add the photo eyepiece.
3. Slide the connecting tube over them and secure it.
4. Screw together the 0.32x camera objective and the adapter piece.
5. Screw on to the adapter piece the camera adapter (T-mount) appropriate to the model of camera.
6. Fit the camera in position.

### Fitting a TV- or cine camera

1. Secure the eyepiece tube in the phototube.
2. Add the photo eyepiece.
3. Mount the  $f=50\text{mm}$  objective with C-mount.
4. Screw the TV- or cine camera into position on top of it.

For the version without photo eyepiece, insert the adapter with C-mount directly into the phototube.

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**Accessory tubes:  
Drawing tube,  
double-iris  
diaphragm,  
discussion tube**

Directions for fitting are on page 24.

**Drawing tube**

- Move the lever to the open circle.
- Position the object and focus on it.
- Move the lever to the white spot.
- Lay the drawing paper under the mirror.
- *Illuminate the drawing paper.* In general, the drawing surface must be illuminated very strongly, but the object only weakly.
- Look into the binocular tube.
- *Balance the illumination brightnesses* so that the object, the drawing surface and the pencil are all seen together.
- Displace the objective so that the drawing surface is focused.

**Accessoires:  
Tube de dessin,  
double diaphragme  
iris, tube de  
discussion**

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


M10 con tubo fotográfico

### Double-iris diaphragm

The lever is used to increase the *depth of field* for spatial objects.

1: Greatest depth of field  
OPEN: Full aperture.

### Discussion tube

Down the discussion tube *two observers* see the same erect, laterally-correct, three-dimensional image.

Please read the separate user manual M2 261.

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
## Accessory tubes: Beam splitters

- 50% observation for TV and dual-station viewing (50% light to each of the two lateral outlets).
- 20% observation for photography and filming (80% light to each of the two lateral outlets).

The double-iris diaphragm affects both the binocular tube and the accessories.

### Fitting

- Fit the beam splitter between the optics carrier and the binocular tube (page 17).
- Attach the accessory to one of the lateral outlets. The accessory is self-centring.

## Second-observer tubes

The monocular and stereo tubes for a second observer can be turned as required.

Use the black rilled sleeve to match the image position to the direction of observation.

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M10 with beam splitter, short tube for second observer, and phototube with Leica camera

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### Phototube

1. Attach the camera adapter to the SLR camera back.
2. Screw the  $f=250\text{mm}$  adapter (large fields) or the  $f=350\text{mm}$  adapter (high image scales) to the camera adapter.
3. Attach the camera and adapter together to the phototube and tighten the clamping screw.
4. Fit the format-indicating graticule for 35mm cameras in one of the eyepieces (page 33).

### Cine/TV tube



1. Screw the C-mount adapter to the cine or TV camera.
2. Attach the camera with adapter to the cine/TV tube and tighten the clamping screw.
3. Slacken the clamping screw to turn the outfit to any of four positions in accordance with requirements. Retighten the clamping screw.

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

M10 with beam splitter and cine/TV attachment

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The stages fit on the incident-light stand with the *spacing ring diam. 100/80mm*. They fit directly into the *80mm diameter port* of the transmitted-light stands and of the stage carriers with magnetic linkage.

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### Gliding stage

The gliding stage can be moved in *all directions* and rotated completely.

- Fit the gliding stage and secure it with the *eccentric lever*.
- Add the glass stage plate or the black/white stage plate, and the stage clips.

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