Cup stage

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Stage carrier with mechanical stage




You can secure the object with the clip and displace it in the x- and y-directions using the rotary knobs (range of movement 76mm x 28mm).

- Fit the stage carrier and immobilize it with the eccentric lever.
- Add the glass insert.
- Using a hollow screw, attach the mechanical stage to the stage carrier.

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Transmitted-light stage

The mirror can be turned to give the most favourable effect in transmitted light.

The reverse side of the mirror produces softer illumination.



- Fit the transmitted-light stage and secure it with the eccentric lever.
- Add the glass stage plate and the stage clips.

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- Connect the 6V/10W low-voltage lamp, complete with lampholder, to the peg, or guide light from another source through the opening and on to the mirror.
- Match the tilt of the mirror to the characteristics of the object.

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


Rotatable polarization stage

- Swing in the polarizer (pull the lever to the left).
- Pull out the compensator.
- Switch on the light.
- Look into the binocular tube and move the lever until the field of view (without the object) is as dark as possible.
- Release the clamping screw and rotate the stage.
- Push in the sensitive tint plate and use the milled ring to turn it.

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
Rotatable Pol. stage with mechanical stage, sensitive-tint plate, and analyzer.

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Assembly

If you have a *sensitive tint plate*:

- Release the screw on the underside of the stage.
- Push the sensitive tint plate fully in.
- Tighten the screw.

If you have a *mechanical stage*:

- Position the mechanical stage so that it can be attached to the rotatable Pol. stage by means of the hollow screw.

- Position the *rotatable Pol. stage* on the transmitted-light stand so that the two red centring screws are aligned symmetrically with respect to the column.
- Tighten the *eccentric lever*.

- Screw the *analyser* to the mount of the objective.

Fit the graticule

1. Unscrew the sleeve from the eyepiece.
2. Push the *graticule* on to the sleeve, ensuring that it sits correctly.
3. Replace the sleeve in the eyepiece.
4. Adjust the *dioptric setting* for the graticule (page 13).

Centring

- Release the clamping screw.
- Turn the stage. The point about which it rotates must coincide with the centre of the crosshair.
- If these points do not coincide, successively displace the red centring screws (Z), experimentally spinning the stage each time, until they do.

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General: Description

Stereomicroscopes impart a spatial image of three-dimensional objects. As in normal observation, each eye observes the object from a different direction. The brain processes the two part-images and fuses them into a single stereoscopic image.

The optical system of the WILD M10 Stereomicroscope consists of a single main objective, behind which are two beam paths which run parallel to one another. The advantages of this design are that each part-image is uniformly sharp from left to right and that the eyes are relaxed and only minimally convergent.

Care

You are working with a high-precision, high-performance instrument, from which you have high expectations.

We therefore guarantee the quality of each of our instruments. Our guarantee covers all faults in materials and manufacture. It does not, however, cover damage resulting from careless or improper handling.

Please handle your valuable optical instrument with respect. It will then fulfil its function with the same precision for decades. Our instruments are famous for that!

If, contrary to expectations, your stereomicroscope develops a fault, please call in a specialist – contact your local Leica agent or the factory in Heerbrugg, Switzerland.

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Dust and dirt will affect the quality of your results

Therefore:

- Put a *dust cover* over the stereomicroscope during breaks in work.
- Use *dust caps* to protect tube openings, tubes without eyepieces, and eyepieces.
- Remove dust with a *pneumatic rubber bulb* or with a *soft brush*.
- Clean objectives and eyepieces with special *optic cleaning cloths* and with pure alcohol.
- When accessories are not being used, *protect them against dust*.

Protect your instruments

- *from damp*, vapours, acids, alkalis, and corrosive substances, and *keep chemicals away* from the instruments.
- *from misuse*.
Do not use non-standard plugs. Never dismantle optical systems or mechanical parts unless instructions for doing so are given in the user manual.
- *from oil and grease*.
Do not grease guide surfaces and mechanical parts.

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What does that mean?

Apochromatically corrected

Chromatic aberrations (colour fringes) are corrected for three points in the spectrum instead of for only two.

Colour temperature

This determines the spectral distribution of the light from a particular source.

Depth of field

The height difference within which features of an object appear acceptably sharp when viewed with an optical instrument.

Dioptric correction

Compensation of long- or short-sightedness by adjusting the eyepieces of the eyepieces.

Distortion-free

The elimination of aberrations which would make a square, flat object look like a cushion or a barrel.

Eyepoint

Rays from all points in the field of view come together at this point. This is where the eye must be positioned.

Field of view

The amount of the object that can be seen with a particular optics combination.

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Focusing

Adjustment of the working distance so that the object is sharply imaged.

Imaging beam path

The path of the image-forming rays, from the object, through the instrument, right to the final image.

Interpupillary distance

The distance between the centres of the pupils of the two eyes (also known as interocular distance).

Parallax-free

Graticule (reticle) and object are sharply imaged in the same plane; this ensures sharp photographs.

Parfocal

The object can be observed from the lowest magnification to the highest without having to refocus.

Plano objective

An objective with a completely-flat field of view.

Polarization

An optical technique for studying objects which are double-refracting (birefringent).

Working distance

The distance from a reference point in the object plane to the front mount of the objective.

Zoom

A continuous magnification changer (as distinct from one which can only be adjusted in steps).

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You can use the table on page 40 to determine the total magnification and field diameter for any magnification changer position and for any optics combination.

After you have correctly set the scale on the magnification changer (page 11), you can read off the approximate total magnification for the following optics combinations:

| Planapo obj. | Eyepiece |
|--------------|---------------|
| 1.0× | 10×, 20×, 32× |
| 0.63× | 10×, 20×, 32× |
| 1.6× | 10×, 20× |

With objectives having focal lengths from 100mm to 350mm

Example

Zoom position 8, objective with focal length 150mm, 10×/21 eyepiece

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