WILD M10

W. NUHSBAUM, INC.

OPTICAL SCIENTIFIC INSTRUMENTS


3918 West Main Street

McHenry, Illinois 60050

(815) 385-5200 FAX (815) 385-5236

User Manual • Mode d'emploi
Bedienungsanleitung • Manual de empleo

Leica

WILD LEITZ ERNEST / GLDESIGN

To you, the user

Thank you for your attention. We wish you pleasure in your future work with this fine stereomicroscope. We are very proud of it, because to begin with all we had was the idea of building a perfect stereomicroscope, using all of the facilities available to us today, and the best and most costly types of glass, without regard to the effort required.

This idea evolved into the creation of the most perfect stereomicroscope the world has ever seen, with the finest optical system and the most modern ergonomics, and with incomparable versatility of accessories. Because top stereomicroscopes come from Wild.

You, as a client of Leica, have available to you the combined expertise of Wild Heerbrugg, Leitz Wetzlar, Cambridge Instruments, Reichert and Jung, and Kern Aarau, when it comes to solving problems. This means that, wherever in the world you happen to be living, you can benefit from personal attention and competent advice – after you have made your purchase, as well.

An efficient after-sales service is at your disposal, if you should need it – which is not likely to happen in a hurry, because instruments from Leica are renowned for their smooth functioning even under extreme conditions or after years of use. What is more, their quality has been confirmed internationally by the award of the coveted ISO 9001 certificate from the International Standards Organization.

A stereomicroscope from Leica really is something special.

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
8

Instructions

- 10 Binocular tube:
Set the viewing height
and the interpupillary
distance
- 10 Eyepieces: Locate the
eyepoint and adjust
the eyecups
- 11 Zoom magnification
changer: Set the
scale, change the
magnification, and
engage the steps
- 12 Adjust the dioptric
setting

Assembly

- 14 Basic outfits with
incident- and
transmitted-light
stands
- 16 Swinging-arm and
table-clamp stands
- 18 Universal stand

Illuminators

- 20 6V/20W low-voltage
illuminator
- 22 Coaxial illuminator

Accessories

- 24 Fitting the accessory
tubes
- 24 Filter-slide housing
- 25 Phototube:
Photography, TV,
filming
- 26 Drawing tube
- 27 Double-iris
diaphragm: Adjust
the depth of field
- 27 Discussion tube:
Observe together
- 28 Beam splitter:
Dual-station viewing,
photography, TV,
filming
- 30 Stages
- 32 Polarization

General

- 34 Description
- 34 Care and maintenance
- 36 Terms used in optics
- 38 Calculating total
magnification and
field diameter
- 40 Optical data

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
Controls on the stereomicroscope

In this section the various functions are *briefly* explained.

Detailed explanations are to be found in the sections dealing with the functions concerned.

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Displace the tubes, using both hands.

Adjust the diopter setting

Rotate the diopter ring on each eyepiece.

Set the magnification scale

Displace the scale in accordance with the optics combination.

Change the magnification

Turn the rotary knob to the magnification setting required.

Set the click-stops

Turn the knurled knob so that the magnification changer engages at nine positions.

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

**Viewing height,
interpupillary
distance,
eye contact,
and magnification
changer**

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Set the viewing height

Tilt the 10°–50° binocular tube upwards or downwards until your head position is comfortable.

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
Adjust the interpupillary distance

Look into the eyepieces and use both hands to displace the eyepiece tubes until you can see a single circular field of view.

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
Locate the eyepoint

Move your head slowly towards the eyepieces until you can see the entire field of view without restriction.

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At this position, the eyepoint, light rays emanating from all points in the field of view come together after leaving the eyepiece. The eyepoint is located 22mm above the eyelens of the 10×/21B wide-field eyepiece for spectacle wearers.

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